

IL17 receptor C antibody

Catalog No: #22229

Package Size: #22229 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	IL17 receptor C antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 62 and 471 of Human IL17RC
Conjugates	Unconjugated
Target Name	IL17 receptor C
Accession No.	Swiss-Prot:Q8NAC3Gene ID:84818
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 86kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug whole cell lysate)
A: HeLa S3
7.5% SDS PAGE
Primary antibody diluted at 1: 1000

Background

This gene encodes a single-pass transmembrane protein that shares limited similarity with the interleukin-17 receptor. Multiple alternatively spliced transcript variants encoding different isoforms have been detected for this gene, but the full-length nature of only three have been determined to date.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.