

EXOC7 antibody

Catalog No: #22233

Package Size: #22233 100ul

Description


Product Name	EXOC7 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 442 and 666 of Human EXOC7
Conjugates	Unconjugated
Target Name	EXOC7
Accession No.	Swiss-Prot:Q9UPT5Gene ID:23265
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 78kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: Molt-4
7.5% SDS PAGE
Primary antibody diluted at 1: 1000

Background

EXOC7 is a component of the exocyst, which is an evolutionarily conserved octameric protein complex essential for exocytosis. The exocyst targets secretory vesicles at specific domains of the plasma membrane for cell surface expansion and protein secretion (Zuo et al., 2006 [PubMed])

Note: This product is for in vitro research use only and is not intended for use in humans or animals.