

ZNF277 antibody

Catalog No: #22242

Package Size: #22242 100ul

Description


Product Name	ZNF277 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 99 and 418 of Human ZNF277
Conjugates	Unconjugated
Target Name	ZNF277
Other Names	NRIF4; ZNF277P
Accession No.	Swiss-Prot:Q9NRM2Gene ID:11179
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 53kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: 293T
10% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.