

PRPS2 antibody

Catalog No: #22246

Package Size: #22246 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PRPS2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 253 of PRPS2
Conjugates	Unconjugated
Target Name	PRPS2
Other Names	PRS II
Accession No.	Swiss-Prot:P11908Gene ID:5634
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 35kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:100-1:500

Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)
A: Hep G2
10% SDS PAGE
PRPS2 antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Hep3B xenograft, using PRPS2 antibody at 1: 100 dilution.



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using PRPS2 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.