

TBRG4 antibody

Catalog No: #22248

Package Size: #22248 100ul

Description

Product Name	TBRG4 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 280 of Human TBRG4
Conjugates	Unconjugated
Target Name	TBRG4
Accession No.	Swiss-Prot:Q969Z0Gene ID:9238
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 58kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:50-1:500

Immunofluorescence: 1:100-1:200


Images




Sample(30 ug whole cell lysate)
A: A431
B: Hep G2
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SW480 xenograft, using TBRG4 antibody at 1: 500 dilution.



Costained with Hoechst 33342



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using TBRG4 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.