

SRPX antibody

Catalog No: #22252

Package Size: #22252 100ul

Description

Product Name	SRPX antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 231 and 434 of SRPX
Conjugates	Unconjugated
Target Name	SRPX
Accession No.	Swiss-Prot:P78539Gene ID:8406
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 52kd


Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images



Sample (30 ug of whole cell lysate)
A: HeLa
10% SDS PAGE
Primary antibody diluted at 1: 1000



Costained with Hoechst 33342



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using SRPX antibody at 1: 500 dilution.

Background

May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.