

Glypican 5 antibody

Catalog No: #22255

Package Size: #22255 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Glypican 5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 188 and 520 of Glypican 5
Conjugates	Unconjugated
Target Name	Glypican 5
Accession No.	Swiss-Prot:P78333Gene ID:2262
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 64kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: NT2D1
7.5% SDS PAGE
Glypican 5 antibody diluted at 1: 500

Background

Cell surface heparan sulfate proteoglycans are composed of a membrane-associated protein core substituted with a variable number of heparan sulfate chains. Members of the glypican-related integral membrane proteoglycan family (GRIPS) contain a core protein anchored to the cytoplasmic

membrane via a glycosyl phosphatidylinositol linkage. These proteins may play a role in the control of cell division and growth regulation. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.