

FLJ35630 antibody

Catalog No: #22256

Package Size: #22256 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	FLJ35630 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 401 and 661 of human BBS12
Conjugates	Unconjugated
Target Name	FLJ35630
Accession No.	Swiss-Prot:Q6ZW61Gene ID:166379
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 79kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: Hep G2
7.5% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.