

ZNF169 antibody

Catalog No: #22257

Package Size: #22257 100ul

Description


Product Name	ZNF169 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 247 and 532 of Human ZNF169
Conjugates	Unconjugated
Target Name	ZNF169
Accession No.	Swiss-Prot:Q14929Gene ID:169841
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 69kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug whole cell lysate)
A: HeLa S3
7.5% SDS PAGE
Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.