

MGC40579 antibody

Catalog No: #22264

Package Size: #22264 100ul

Description

Product Name	MGC40579 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 297 of Human GK5
Conjugates	Unconjugated
Target Name	MGC40579
Accession No.	Swiss-Prot:Q6ZS86Gene ID:256356
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 34kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Images




Sample (30 ug of whole cell lysate)

A: A431

10% SDS PAGE

Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded CL1-0 xenograft, using GK5 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.