

PODXL2 antibody

Catalog No: #22268

Package Size: #22268 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PODXL2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 273 and 567 of PODXL2
Conjugates	Unconjugated
Target Name	PODXL2
Accession No.	Swiss-Prot:Q9NZ53Gene ID:50512
Concentration	0.8mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 65kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: Hela
7.5% SDS PAGE
PODXL2 antibody diluted at 1: 1000

Background

Belongs to the podocalyxin family. Homodimer; disulfide-linked. Binds to SELL/L-selectin. 2 isoforms of the human protein are produced by alternative splicing.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.