

Edc3 antibody

Catalog No: #22272

Package Size: #22272 100ul

Description

Product Name	Edc3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 127 and 346 of Edc3
Conjugates	Unconjugated
Target Name	Edc3
Accession No.	Swiss-Prot:Q96F86Gene ID:80153
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 56kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: A431
7.5% SDS PAGE
Primary antibody diluted at 1: 3000



Immunohistochemical analysis of paraffin-embedded MCF7 xenograft, using Edc3 antibody at 1: 250 dilution.

Background

EDC3 is associated with an mRNA-decapping complex required for removal of the 5-prime cap from mRNA prior to its degradation from the 5-prime end (Fenger-Gron et al., 2005 [PubMed 16364915]).[supplied by OMIM]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.