

Alpha-catulin antibody

Catalog No: #22278

Package Size: #22278 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha-catulin antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 292 and 533 of Human CTNNAL1
Conjugates	Unconjugated
Target Name	Alpha-catulin
Accession No.	Swiss-Prot:Q9UBT7Gene ID:8727
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 82kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images




Sample (30 ug of whole cell lysate)

A: HeLa

B: Hep G2

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SW480 xenograft, using Alpha-catulin antibody at 1: 100 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.