

PABPC3 antibody

Catalog No: #22285

Package Size: #22285 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PABPC3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 107 and 293 of PABPC3
Conjugates	Unconjugated
Target Name	PABPC3
Accession No.	Swiss-Prot:Q9H361Gene ID:5042
Concentration	0.8mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 70kd

Immunohistochemistry: 1:100-1:250

Images




Sample (30 ug of whole cell lysate)

A: JurKat

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SCM1 xenograft, using PABPC3 antibody at 1: 100 dilution.

Background

Messenger RNA stability and translation initiation are extensively under the control of poly(A)-binding proteins (PABP). See PABPC1 (MIM 604679) for background information.[supplied by OMIM]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.