

OMG antibody

Catalog No: #22290

Package Size: #22290 100ul

Description


Product Name	OMG antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 249 and 301 (P23515) of OMG
Conjugates	Unconjugated
Target Name	OMG
Accession No.	Swiss-Prot:P23515Gene ID:4974
Concentration	1mg/ml
Formulation	Supplied in 1XPBS, 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 50kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: IMR32
10% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.