

ARHGEF5 antibody

Catalog No: #22309

Package Size: #22309 100ul

Description


Product Name	ARHGEF5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Description	A synthesized peptide derived from human ARHGEF5
Conjugates	Unconjugated
Target Name	ARHGEF5
Accession No.	NCBI Gene ID: 7984NCBI mRNA#: BC010046NCBI Protein#: AAH10046
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 60kd

Western blotting: 1:500-1:1000

Images



Western blot analysis of extracts from RAW264.7 cells, using ARHGEF5 antibody.

Background

Rho GTPases play a fundamental role in numerous cellular processes initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein may be involved in the control of cytoskeletal organization. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.