

ARF5 antibody

Catalog No: #22315

Package Size: #22315 100ul

Description


Product Name	ARF5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 49 and 113 of Human ARF5
Conjugates	Unconjugated
Target Name	ARF5
Accession No.	Swiss-Prot:P84085Gene ID:381
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 21kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: MOLT4
B: Raji
15% SDS PAGE
Primary antibody diluted at 1: 1500

Background

ADP-ribosylation factor 5 (ARF5) is a member of the human ARF gene family. These genes encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking and as activators of phospholipase D. The gene

products include 6 ARF proteins and 11 ARF-like proteins and constitute 1 family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2, and ARF3), class II (ARF4 and ARF5) and class III (ARF6). The members of each class share a common gene organization. The ARF5 gene spans approximately 3.2kb of genomic DNA and contains six exons and five introns. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.