

TENC1 antibody

Catalog No: #22371

Package Size: #22371 100ul

Description


Product Name	TENC1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1098 and 1419 of TENC1
Conjugates	Unconjugated
Target Name	TENC1
Other Names	C1-TEN; C1TEN; DKFZp686D13244; FLJ16320; KIAA1075; TNS2
Accession No.	Swiss-Prot:Q63HR2Gene ID:23371
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 154kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: HeLa
5% SDS PAGE
Primary antibody diluted at 1: 1000

Background

The protein encoded by this gene belongs to the tensin family. Tensin is a focal adhesion molecule that binds to actin filaments and participates in

signaling pathways. This protein plays a role in regulating cell migration. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.