

ACAD8 antibody

Catalog No: #22378

Package Size: #22378 100ul

Description


Product Name	ACAD8 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 146 and 411 of Human ACAD8
Conjugates	Unconjugated
Target Name	ACAD8
Other Names	ACAD-8; FLJ22590
Accession No.	Swiss-Prot:Q9UKU7Gene ID:27034
Concentration	0.3mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 45kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: 293T
10% SDS PAGE
Primary antibody diluted at 1: 500

Background

Isobutyryl-CoA dehydrogenase (EC 1.3.99), encoded by the ACAD8 gene, catalyzes the third step of the degradation of the branched chain amino

Note: This product is for in vitro research use only and is not intended for use in humans or animals.