

HPS6 antibody

Catalog No: #22400

Package Size: #22400 100ul

Description


Product Name	HPS6 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 405 and 456 (Q86YV9) of HPS6
Conjugates	Unconjugated
Target Name	HPS6
Accession No.	Swiss-Prot:Q86YV9Gene ID:79803
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 83kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: A431

7.5% SDS PAGE

Primary antibody diluted at 1: 1000

Background

This intronless gene encodes a protein that may play a role in organelle biogenesis associated with melanosomes, platelet dense granules, and lysosomes. This protein interacts with Hermansky-Pudlak syndrome 5 protein. Mutations in this gene are associated with Hermansky-Pudlak

syndrome type 6. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.