

## TCP-1 theta antibody

Catalog No: #22412

Package Size: #22412 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                        |
|-----------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name          | TCP-1 theta antibody                                                                                                   |
| Host Species          | Rabbit                                                                                                                 |
| Clonality             | Polyclonal                                                                                                             |
| Purification          | Purified by antigen-affinity chromatography.                                                                           |
| Applications          | WB IHC                                                                                                                 |
| Species Reactivity    | Hu                                                                                                                     |
| Immunogen Type        | Recombinant protein                                                                                                    |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 111 and 337 of Human CCT8 |
| Conjugates            | Unconjugated                                                                                                           |
| Target Name           | TCP-1 theta                                                                                                            |
| Accession No.         | Swiss-Prot:P50990Gene ID:10694                                                                                         |
| Concentration         | 1mg/ml                                                                                                                 |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.         |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                              |


## Application Details

Predicted MW: 60kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

## Images



Sample (30 ug of whole cell lysate)


A: A431

B: H1299

C: HeLa

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded OVCAR3 xenograft, using TCP-1 theta antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.