

GPCRC6A antibody

Catalog No: #22439

Package Size: #22439 100ul

Description


Product Name	GPCRC6A antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 862 and 926 of Human GPCR6A
Conjugates	Unconjugated
Target Name	GPCRC6A
Other Names	GPCR; bA86F4.3
Accession No.	Swiss-Prot:Q5T6X5Gene ID:222545
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 105kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: Raji
7.5% SDS PAGE
Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.