

RAD51AP1 antibody

Catalog No: #22444

Package Size: #22444 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	RAD51AP1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 273 and 335 of Human RAD51AP1
Conjugates	Unconjugated
Target Name	RAD51AP1
Accession No.	Swiss-Prot:Q96B01Gene ID:10635
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 37kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)

A: Raji

12% SDS PAGE

Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.