

PLEK antibody

Catalog No: #22455

Package Size: #22455 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PLEK antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 287 and 350 of Human PLEK
Conjugates	Unconjugated
Target Name	PLEK
Other Names	P47; FLJ27168
Accession No.	Swiss-Prot:P08567Gene ID:5341
Concentration	0.6mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 40kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)

A: Raji

12% SDS PAGE

Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.