

Profilin 2 antibody

Catalog No: #22457

Package Size: #22457 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Profilin 2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 120 of Human PFN2
Conjugates	Unconjugated
Target Name	Profilin 2
Accession No.	Swiss-Prot:P35080Gene ID:5217
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 15kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug whole cell lysate)
A: A431
B: H1299
12% SDS PAGE
Primary antibody diluted at 1: 500

Background

The protein encoded by this gene is a ubiquitous actin monomer-binding protein belonging to the profilin family. It is thought to regulate actin polymerization in response to extracellular signals. There are two alternatively spliced transcript variants encoding different isoforms described for this

gene. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.