

BLU antibody

Catalog No: #22462

Package Size: #22462 100ul

Description


Product Name	BLU antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 377 and 440 of Human ZMYND10
Conjugates	Unconjugated
Target Name	BLU
Accession No.	Swiss-Prot:O75800Gene ID:51364
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 50kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)

A: 239T

B: A431

C: H1299

12% SDS PAGE

Primary antibody diluted at 1: 500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.