

GIPC1 antibody

Catalog No: #22473

Package Size: #22473 100ul

Description

Product Name	GIPC1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 273 and 333 of Human GIPC1
Conjugates	Unconjugated
Target Name	GIPC1
Accession No.	Swiss-Prot:O14908Gene ID:10755
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 36kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: A431
B: HeLa
10% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SCC4 xenograft, using GIPC1 antibody at 1: 100 dilution.

Background

GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking (Lee et al., 2008 [PubMed 18775991]).[supplied by OMIM]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.