

SAMHD1 antibody

Catalog No: #22479

Package Size: #22479 100ul

Description

Product Name	SAMHD1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 122 and 421 of SAMHD1
Conjugates	Unconjugated
Target Name	SAMHD1
Other Names	DCIP
Accession No.	Swiss-Prot:Q9Y3Z3Gene ID:25939
Concentration	0.6mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 72kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: Raji
7.5% SDS PAGE
SAMHD1 antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded OVCA, using SAMHD1 antibody at 1: 100 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.