

PGAM1 antibody

Catalog No: #22490

Package Size: #22490 100ul


Description

Product Name	PGAM1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 16 and 209 of PGAM1
Conjugates	Unconjugated
Target Name	PGAM1
Accession No.	Swiss-Prot:P18669Gene ID:5223
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 29kd

Images



Sample (30 ug of whole cell lysate)

A: 293T

B: A431

C: JurKat

D: Raji

12% SDS PAGE

Primary antibody diluted at 1: 1000

Background

Interconversion of 3- and 2-phosphoglycerate with 2,3-bisphosphoglycerate as the primer of the reaction. Can also catalyze the reaction of EC 5.4.2.4 (synthase) and EC 3.1.3.13 (phosphatase), but with a reduced activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.