

BZW2 antibody

Catalog No: #22515

Package Size: #22515 100ul

Description

Product Name	BZW2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 208 of BZW2
Conjugates	Unconjugated
Target Name	BZW2
Accession No.	Swiss-Prot:Q9Y6E2Gene ID:28969
Concentration	0.6mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 48kd


Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200


Images



Sample(30 ug whole cell lysate)
A: A431
B: HeLa S3
C: MOLT4
10% SDS PAGE
Primary antibody diluted at 1: 1000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using BZW2 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.