

PIK3R3 antibody

Catalog No: #23086

Package Size: #23086 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PIK3R3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 192 and 437 of Human PIK3R3
Conjugates	Unconjugated
Target Name	PIK3R3
Accession No.	Swiss-Prot:Q92569Gene ID:8503
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 54kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Images



Sample (30ug whole cell lysate)
A: A431
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SAS xenograft, using PIK3R3 antibody at 1: 500 dilution.

Published Papers

Yan Xia Yu;Hai Jian Wu;Bing Xu Tan;Chen Qiu;Hui Zhong Liu et al., CSF-1R regulates non-small cell lung cancer cells dissemination through Wnt3a signaling., , (2017)

[PMID:29218239](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.