

RILP Antibody


Catalog No: #25472

Package Size: #25472 100ul

Description

Product Name	RILP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB
Species Reactivity	Hu Ms
Immunogen Type	Peptide
Immunogen Description	RILP antibody was raised against an 18 amino acid peptide near the center of human RILP.
Conjugates	Unconjugated
Target Name	RILP
Other Names	Rab-interacting lysosomal protein, PP10141
Accession No.	Swiss-Prot:Q96NA2Gene ID:83547
Concentration	1mg/ml
Formulation	PBS containing 0.02% sodium azide.
Storage	Store at 4°C for three months and -20°C, stable for up to one year.

Images



Western blot analysis of RILP in HeLa cell lysate with RILP antibody at 1 ug/mL in (A) the absence and (B) the presence of blocking peptide.

Published Papers

Hai Qiu;Manting Zhang;Minchun Li;Chuanxiang Chen;Huijun Wang;Xia Yue et al., Methamphetamine and Methamphetamine-Induced Neuronal Exosomes Modulate the Activity of Rab7a via PTEN to Exert an Influence on the Disordered Autophagic Flux Induced in Neurons., (2025)
[PMID:40141286](https://pubmed.ncbi.nlm.nih.gov/40141286/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.