

TBL1Y Monoclonal Antibody

Catalog No: #27006

Package Size: #27006-1 50ul #27006-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	TBL1Y Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms P Rt D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant full length protein (522aa) of human TBL1Y.
Conjugates	Unconjugated
Target Name	TBL1Y
Other Names	F-box-like/WD repeat-containing protein TBL1Y ☐ transducin beta-like 1 ☐ transducin beta-like protein 1Y ☐ transducin-beta-like protein 1, Y-linked
Accession No.	UniProt:Q9BQ87 NCBI Gene ID:90665
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 60kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected TBL1Y on 40ug of cell lysates



1:150 dilution from a previous lot detected TBL1Y in HeLa

Note: This product is for in vitro research use only and is not intended for use in humans or animals.