

CBX5 Monoclonal Antibody

Catalog No: #27012

Package Size: #27012-1 50ul #27012-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CBX5 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Ms P Rt D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant full length protein (191aa) of human CBX5 .
Conjugates	Unconjugated
Target Name	CBX5
Other Names	chromobox protein homolog 5 □HP1-ALPHA □HP1Hs alpha □antigen p25 □HP1 alpha homolog □heterochromatin protein 1-alpha □heterochromatin protein 1 homolog alpha □chromobox homolog 5 (HP1 alpha homolog, Drosophila)
Accession No.	UniProt:P45973 NCBI Gene ID:23468
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C


Application Details

Predicted MW: 26kd
Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected CBX5 on 40ug of cell lysates



1:100 dilution from a previous lot detected CBX5 in HeLa

Note: This product is for in vitro research use only and is not intended for use in humans or animals.