

P53 Monoclonal Antibody

Catalog No: #27021

Package Size: #27021-1 50ul #27021-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	P53 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu D Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment of human TP53 (C-terminus)
Conjugates	Unconjugated
Target Name	P53
Other Names	cellular tumor antigen p53 □ antigen NY-CO-13 □ phosphoprotein p53 □ p53 tumor suppressor □ transformation-related protein 53
Accession No.	UniProt:P04637 NCBI Gene ID:7157
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 53kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected TP53 on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.