

TRIM28 Monoclonal Antibody

Catalog No: #27026

Package Size: #27026-1 50ul #27026-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TRIM28 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC ICC
Species Reactivity	Hu D Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 548aa-825aa of human TRIM28(C-terminus).
Conjugates	Unconjugated
Target Name	TRIM28
Other Names	transcription intermediary factor 1-beta □KAP-1 □KRIP-1 □TIF1-beta □RING finger protein 96 □KRAB-associated protein 1 □nuclear corepressor KAP-1 □KRAB-interacting protein 1 □E3 SUMO-protein ligase TRIM28 □tripartite motif-conta
Accession No.	UniProt:Q13263 NCBI Gene ID:10155
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C


Application Details

Predicted MW: 100kd
Western blotting: 1:1000


Images



1:4000 dilution of this antibody detected TRIM28 on 40ug of cell lysates



HeLa cells using anti-TRIM28 (C-terminus)antibody diluted 1:200



Paraffin-embedded human breast cancer using anti-TRIM28 (C-terminus) diluted1:50 -1:100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.