

Albumin Monoclonal Antibody

Catalog No: #27044

Package Size: #27044-1 50ul #27044-2 100ul

Description


Product Name	Albumin Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human Albumin
Conjugates	Unconjugated
Target Name	Albumin
Other Names	serum albumin □ albumin (32 AA) □ albumin (AA 34) □ growth-inhibiting protein 20 □ cell growth inhibiting protein 42
Accession No.	UniProt:P02768 NCBI Gene ID:213
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 68kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected Albumin on 40ug of cell lysates



Paraffin--embedded human liver cancer using anti- ALBdiluted
1:50-1:100

Published Papers

el at., Elevated bioaccumulation of PFAAs in *Oryzias melastigma* following the increase of salinity is associated with the up-regulated expression of PFAA-binding proteins. In *Sci Total Environ* on 2020 Apr 3; by Avell n-Llaguno RD, Liu X, et al.. PMID:32298882 , , (2020)

[PMID:32298882](#)

el at., Strontium ions protect hearts against myocardial ischemia/reperfusion injury. In *Sci Adv* on 2021 Jan 15 by Min Xing, Yun Jiang, et al.. PMID:33523909, , (2021)

[PMID:33523909](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.