

DsbA Monoclonal Antibody

Catalog No: #27049

Package Size: #27049-1 50ul #27049-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	DsbA Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	E.coli
Immunogen Type	Recombinant protein
Immunogen Description	Full Length (Escherichia coli) Disulfide Oxidoreductase, isolated from E. coli.
Conjugates	Unconjugated
Target Name	DsbA
Other Names	hypothetical protein
Accession No.	UniProt:C5WBA2 NCBI Gene ID:8114936
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

Application Details

Predicted MW: 24kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected DsbA on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.