

Ikaros(C-term) Monoclonal Antibody

Catalog No: #27062

Package Size: #27062-1 50ul #27062-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Ikaros(C-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms D
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 249-415aa of human IKAROS(C-terminus).
Conjugates	Unconjugated
Target Name	Ikaros
Other Names	DNA-binding protein Ikaros □CLL-associated antigen KW-6 □Ikaros (zinc finger protein) □ikaros family zinc finger protein 1 □lymphoid transcription factor LyF-1 □zinc finger protein, subfamily 1A, 1 (Ikaros
Accession No.	UniProt:Q13422 NCBI Gene ID:10320
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details

Predicted MW: 58kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected IKAROS on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.