

## AMACR(C-term) Monoclonal Antibody

Catalog No: #27068

Package Size: #27068-1 50ul #27068-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                       |                                                                                                                             |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name          | AMACR(C-term) Monoclonal Antibody                                                                                           |
| Host Species          | Mouse                                                                                                                       |
| Clonality             | Monoclonal                                                                                                                  |
| Isotype               | IgG2b                                                                                                                       |
| Purification          | Purified by antigen-affinity chromatography.                                                                                |
| Applications          | WB ICC                                                                                                                      |
| Species Reactivity    | Hu Ms Rt                                                                                                                    |
| Immunogen Type        | Recombinant protein                                                                                                         |
| Immunogen Description | Recombinant fragment, corresponding to amino acids 138-319aa of human AMACR (C-terminus).                                   |
| Conjugates            | Unconjugated                                                                                                                |
| Target Name           | AMACR                                                                                                                       |
| Other Names           | alpha-methylacyl-CoA racemase □ 2-methylacyl-CoA racemase                                                                   |
| Accession No.         | UniProt:Q9UHK6<br>NCBI Gene ID:23600                                                                                        |
| Formulation           | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol |
| Storage               | store at -20°C                                                                                                              |

## Application Details


Predicted MW: 42kd

Western blotting: 1:1000

## Images



1:1000 dilution of this antibody detected AMACR on 40ug of cell lysates



HeLa cells using anti- AMACR (C-terminus) antibody diluted 1:150

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.