

CDK4(N-term) Monoclonal Antibody

Catalog No: #27070

Package Size: #27070-1 50ul #27070-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CDK4(N-term) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Purification	Purified by antigen-affinity chromatography.
Applications	WB ICC
Species Reactivity	Hu Rb P D Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids 2-151aa of human CDK4(N-terminus).
Conjugates	Unconjugated
Target Name	CDK4
Other Names	cyclin-dependent kinase 4 □ cell division protein kinase 4
Accession No.	UniProt:P11802 NCBI Gene ID:1019
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Storage	store at -20°C

Application Details


Predicted MW: 34kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected CDK4 on 40ug of cell lysates



HeLa cells using anti- CDK4(N-terminus)antibody diluted
1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.