

Hexokinase-2(N-term) Monoclonal Antibody

Catalog No: #27074

Package Size: #27074-1 50ul #27074-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|-----------------------|--|
| Product Name | Hexokinase-2(N-term) Monoclonal Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Purification | Purified by antigen-affinity chromatography. |
| Applications | WB FC |
| Species Reactivity | Hu P Ms Rt |
| Immunogen Type | Recombinant protein |
| Immunogen Description | Recombinant fragment, corresponding to amino acids 270-435aa of human HK2(N-terminus). |
| Conjugates | Unconjugated |
| Target Name | Hexokinase-2 |
| Other Names | hexokinase-2 <input type="checkbox"/> HK II <input type="checkbox"/> hexokinase type II <input type="checkbox"/> hexokinase-2, muscle <input type="checkbox"/> muscle form hexokinase <input type="checkbox"/> mi-2 autoantigen 240 kDa protein <input type="checkbox"/> zinc-finger helicase (Snf2-like) |
| Accession No. | UniProt:P52789 NCBI Gene ID:3099 |
| Formulation | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol |
| Storage | store at -20°C |

Application Details

Predicted MW: 102kd


Western blotting: 1:1000

Images



1:3000 dilution of this antibody detected HK2 on 40ug of cell lysates

K562 cells stained with HK2 (red, 1/100 dilution), followed by FITC-conjugated goat anti-mouse IgG. Blue line histogram represents the isotype control, normal mouse IgG



Note: This product is for in vitro research use only and is not intended for use in humans or animals.