

ARRB1 Monoclonal Antibody

Catalog No: #27077

Package Size: #27077-1 50ul #27077-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	ARRB1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant fragment, corresponding to amino acids of human ARRB1
Conjugates	Unconjugated
Target Name	ARRB1
Other Names	arrestin-beta-1
Accession No.	UniProt:P49407 NCBI Gene ID:408
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Storage	store at -20°C

Application Details

Predicted MW: 47kd

Western blotting: 1:1000

Images



1:1000 dilution of this antibody detected ARRB1 on 40ug of cell lysates

Note: This product is for in vitro research use only and is not intended for use in humans or animals.