

Brgl Monoclonal Antibody

Catalog No: #27090

Package Size: #27090-1 50ul #27090-2 100ul

Description


Product Name	Brgl Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4E11-E11-G11
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human SMARCA4 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	Brgl
Other Names	BRG1, SNF2, SWI2, RTPS2, BAF190, SNF2L4, SNF2LB, hSNF2b, SNF2-BETA
Accession No.	Swiss-Prot:P51532Gene ID:6597
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Application Details

Predicted MW: 235kd

Western blotting: 1:1000

Images



Western blot detection of SMARCA4 in K562&HeLaNE cell lysates using SMARCA4 antibody (1:1000 diluted).
Predicted band size: 235KDa
Observed band size: 235KDa.

Background

Transcriptional coactivator cooperating with nuclear hormone receptors to potentiate transcriptional activation. Component of the CREST-BRG1 complex, a multiprotein complex that regulates promoter activation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.