

SAFB-1 Monoclonal Antibody

Catalog No: #27104

Package Size: #27104-1 50ul #27104-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SAFB-1 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2E8-E2-G6
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human SAFB-1 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	SAFB-1
Other Names	HAP, HET, SAFB1, SAF-B1
Accession No.	Swiss-Prot:Q9BU40Gene ID:91851
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol


Application Details

Predicted MW: 130kd
Western blotting: 1:1000

Images



Western blot detection of SAFB-1 in HeLaNE,A431,MDA-MB-468,Breast cancer,ZR75-1&MCF7 cell lysates using SAFB-1 antibody (1:4000 diluted).
Predicted band size: 130kDa
Observed band size: 130kDa



Immunocytochemistry stain of HeLa using SAFB-1 antibody (1:300).

Background

Binds to scaffold/matrix attachment region (S/MAR) DNA and forms a molecular assembly point to allow the formation of a 'transcriptosomal' complex (consisting of SR proteins and RNA polymerase II) coupling transcription and RNA processing. By similarity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.