

NUP98 Monoclonal Antibody

Catalog No: #27111

Package Size: #27111-1 50ul #27111-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NUP98 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3B8-D7-H10
Isotype	IgG1
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu D Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human NUP98 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	NUP98
Other Names	ADIR2, NUP96, NUP196
Accession No.	Swiss-Prot:P31939Gene ID:471
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol


Application Details

Predicted MW: 98kd
Western blotting: 1:1000

Images



Western blot detection of NUP98 in HeLa



Immunoprecipitation analysis of HeLa cell lysates using NUP98 antibody

Background

Nup98 and Nup96 play a role in the bidirectional transport across the nucleoporin complex (NPC). The repeat domain in Nup98 has a direct role in the transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.