

TXNRD2 Monoclonal Antibody

Catalog No: #27123

Package Size: #27123-1 50ul #27123-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	TXNRD2 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3F2-E12-F10
Isotype	IgG1
Applications	ELISA WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human TXNRD2 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	TXNRD2
Other Names	TR, TR3, SELZ, TRXR2, TR-BETA
Accession No.	Swiss-Prot:Q9UH65Gene ID:23075
Formulation	ascites

Application Details

Predicted MW: 56kd

Western blotting: 1:1000

Images



Western blot detection of TXNRD2 in A431,U251

Background

Maintains thioredoxin in a reduced state. Implicated in the defenses against oxidative stress. May play a role in redox-regulated cell signaling.

Maintains thioredoxin in a reduced state. Implicated in the defenses against oxidative stress.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.