

SWAP70 Monoclonal Antibody

Catalog No: #27126

Package Size: #27126-1 50ul #27126-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	SWAP70 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2E9-E8-E9
Isotype	IgG2b
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human SWAP70 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	SWAP70
Other Names	HSPC321, SWAP-70, FLJ39540, KIAA0640
Accession No.	Swiss-Prot:Q86TG7 Gene ID:23089
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Application Details

Predicted MW: 69kd

Western blotting: 1:1000

Images



Western blot detection of SWAP70 in Rat liver

Background

Phosphatidylinositol 3,4,5-trisphosphate-dependent guanine nucleotide exchange factor (GEF) which, independently of RAS, transduces signals from tyrosine kinase receptors to RAC. It also mediates signaling of membrane ruffling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.