

CRK Monoclonal Antibody

Catalog No: #27127

Package Size: #27127-1 50ul #27127-2 100ul

Description


Product Name	CRK Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3H7-E5-H8
Isotype	IgG2b
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Hu D Rb Ms
Immunogen Type	Recombinant protein
Immunogen Description	Purified recombinant human CRK protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	CRK
Other Names	p38, CRKII
Accession No.	Swiss-Prot:P05362Gene ID:3383
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Application Details

Predicted MW: 34kd

Western blotting: 1:1000

Images



Western blot detection of CRK in LNCAP,A431

Background

The Crk-I and Crk-II forms differ in their biological activities. Crk-II has less transforming activity than Crk-I. Crk-II mediates attachment-induced MAPK8 activation, membrane ruffling and cell motility in a Rac-dependent manner.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.