

MPHOSPH6 Polyclonal Antibody

Catalog No: #29416

Package Size: #29416-1 50ul #29416-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MPHOSPH6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MPHOSPH6 (NP_005783.2).
Conjugates	Unconjugated
Other Names	MPHOSPH6; MPP; MPP-6; MPP6; M-phase phosphoprotein 6
Accession No.	Swiss-Prot#:Q99547NCBI Gene ID:10200
Calculated MW	22kDa,23kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MPHOSPH6 at 1:1000 dilution.

Background

Published Papers

el at., 5-(3,4-Difluorophenyl)-3-(6-methylpyridin-3-yl)-1,2,4-oxadiazole (DDO-7263), a novel Nrf2 activator targeting brain tissue, protects against MPTP-Induced subacute Parkinson's disease in mice by inhibiting the NLRP3 inflammasome and protects PC12 cells against oxidative stress. In

Note: This product is for in vitro research use only and is not intended for use in humans or animals.